

**Anti-Human/Mouse/Rat ITGB6 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12358**

---

3th Edition

**Description:** Integrin beta-6 is a protein that in humans is encoded by the ITGB6 gene. Integrin, beta 6 has been shown to interact with FHL2. Integrin alpha-V/beta-6 is a receptor for fibronectin and cytotactin. It recognizes the sequence R-G-D in its ligands. Internalisation of integrin alpha-V/beta-6 via clathrin-mediated endocytosis promotes carcinoma cell invasion.

**Antigen:** Recombinant protein of human ITGB6

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**