

Anti-Human/Mouse CD40LG Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12403

3th Edition

Description:The protein encoded by this gene is expressed on the surface of T cells. It regulates B cell function by engaging CD40 on the B cell surface. A defect in this gene results in an inability to undergo immunoglobulin class switch and is associated with hyper-IgM syndrome

Antigen:Recombinant protein of human CD40LG

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Specifically expressed on activated CD4+ T-lymphocytes.