

**Anti-Human/Mouse/Rat MPZ Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12476**

---

3th Edition

**Description:** This gene encodes a major structural protein of peripheral myelin. Mutations in this gene result in the autosomal dominant form of Charcot-Marie-Tooth disease type 1 and other polyneuropathies. The myelin sheath is a multi-layered membrane, unique to the nervous system, that functions as an insulator to greatly increase the velocity of axonal impulse conduction. The P0 glycoprotein, absent in the central nervous system,[1] is a major component of the myelin sheath in peripheral nerves. It comprises a large extracellular N-terminal domain, a single transmembrane (TM) region, and a smaller positively charged intracellular domain.

**Antigen:** Recombinant protein of human MPZ

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**