

Anti-Human/Mouse/Rat PPP3CA Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12477

3th Edition

Description:Protein phosphatase 3, catalytic subunit, alpha isozyme is a protein that in humans is encoded by thePPP3CA gene. Calcium-dependent, calmodulin-stimulated protein phosphatase. Many of the substrates contain a PxlxIT motif. This subunit may have a role in the calmodulin activation of calcineurin. Dephosphorylates DNM1L, HSPB1 and SSH1.

Antigen:Recombinant protein of human PPP3CA

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: