

Anti-Human/Mouse/Rat HIPK4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12734

3th Edition

Description: HIPK4 (Homeodomain-interacting protein kinase 4) is a 616 amino acid protein that localizes to the cytoplasm and exists as two isoforms that are produced by alternative splicing events. Like other members of the HIPK family, HIPK4 functions as a protein kinase that catalyzes the ATP-dependent phosphorylation of target proteins and is thought to act as a corepressor of target transcription factors.

Antigen: Recombinant protein of human HIPK4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: