

**CNTF(Ciliary neurotrophic factor), 1-200aa, Human, His-tagged, Recombinant, E.coli**

**Cat.NO.: TP01670**

3th Edition

**Synonyms:**Ciliary neurotrophic factor

**Description:**CNTF (Ciliary neurotrophic factor) is a polypeptide hormone whose actions appear to be restricted to the nervous system where it promotes neurotransmitter synthesis and neurite outgrowth in certain neuronal populations. This protein is a potent survival factor for neurons and oligodendrocytes and may be relevant in reducing tissue destruction during inflammatory attacks.

**Form:**Liquid. In 20 mM Tris-HCl buffer (pH 8.5) containing 1 mM DTT, 30% Glycerol, 0.2M NaCl

**Molecular Weight:**25kDa (220aa), confirmed by MALDI-TOF.

**Sequences:**

MGSSHHHHHHSSGLVPRGSHMAFTEHSPLTPHRRDLCRSRWLARKIRSDLTALTESYVKHQGLNKNINLDSADG  
MPVASTDQWSELTEAERLQENLQAYRTFHVLLARLLEDQQVHFTPTGDFHQAIHTLLLQVAAFAYQIEELMILLEY  
KIPRNEADGMPINVGDGGLFEKKLWGLKVLQELSQWTVRSIHDLRFISSHQTGIPARGSHYIANNKKM

**Purity:**> 95% by HPLC

**Concentration:**1 mg/ml (determined by Bradford assay)

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.