

FBP1, 1-338aa, Human, His-tag, E.coli

Cat.NO.: TP02096

3th Edition

Synonyms:FBPase1, FBP, FBP1.

Description:Fructose-1,6-bisphosphatase 1(FBP1) is a gluconeogenesis regulatory enzyme, catalyzes the hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate and inorganic phosphate. Fructose-1,6-diphosphatase deficiency is associated with hypoglycemia and metabolic acidosis. Recombinant FBP1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

Form:Liquid. In 20 mM Tris-HCl buffer (pH8.0) containing 1mM DTT 10% glycerol

Molecular Weight:39.0 kDa (358aa)

Sequences:

MGSSHHHHHSSGLVPRGSHMADQAPFDTDVNTLTRFVMEEGRKARGTGELTQLLNSLCTAVKAISSAVRKAGIA
HLYGIAGSTNVTGDQVKKLDVLSNDLVMNMLKSSFATCVLVSEEDKHAIIVEPEKRGKYVVCFDPLDGSSNIDCLVS
VGTIFGIYRKKSTDEPSEKDALQPGRNLVAAGYALYGSATMLVLAMDCGVNCFMLDPAIGEFILVDKDVKIKKKGKIY
SLNEGYARDFDPAVTEYIQRKKFPPDNSAPYGARYVGSMAVADVHRTLVIYGGIFLYPANKKSPNGKLRLLYECPMA
YVMEKAGGMATTGKEAVLDVIPTDIHQRAPVILGSPDDVLEFLKVYEKHSQA

Purity:> 95% by HPLC

Concentration:1 mg/ml (determined by Bradford assay)

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.