

**Recombinant Mouse REG3B / Pap1 Protein (His tag)**

**Cat.NO.: TP07088**

3th Edition

**Synonyms:**HIP;Pap;PAP1;REG-III

**Description:**

**Form:**PBS

**Molecular Weight:**18.1 kDa

**Sequences:**Met1-Gly175

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.