

**Recombinant Human CD3E & CD3G Heterodimer Protein**

**Cat.NO.: TP07097**

3th Edition

**Synonyms:**CD3E & CD3G Heterodimer;IMD18;T3E;TCRE

**Description:**

**Form:**PBS

**Molecular Weight:**77.2 (40.1+37.1) kDa

**Sequences:**Met 1-Asp126

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.