

Recombinant Human CD3D & CD3E Heterodimer Protein

Cat.NO.: TP07098

3th Edition

Synonyms:CD3-DELTA;IMD19;T3D

Description:

Form:PBS

Molecular Weight:77.5 (37.5+40? kDa

Sequences:Met 1-Ala 105

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.