

Recombinant Human Gastric intrinsic factor / GIF Protein (His tag)

Cat.NO.: TP07133

3th Edition

Synonyms:IF;IFMH;INF;TCN3

Description:Gastric intrinsic factor, also known as GIF, belongs to the of the cobalamin transport protein family. It is a glycoprotein produced by the parietal cells of the stomach. Gastric intrinsic factor plays a key role in the absorption of vitamin B12 on in the small intestine. Vitamin B12 bounds to haptocorrin after entry into the stomach. The resulting complex enters the duodenum, where pancreatic enzymes digest haptocorrin. In the less acidic environment of the small intestine, B12 can then bind to gastric intrinsic factor. This new complex travels to the ileum, where special epithelial cells endocytose them. Inside the cell, B12 dissociates once again and binds to another protein, transcobalamin II. The new complex can exit the epithelial cells to enter the liver.

Form:PBS

Molecular Weight:44.8 kDa

Sequences:Met 1-Tyr417

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.