

**Recombinant Human SerpinB3 Protein (His Tag)**

**Cat.NO.: TP07321**

3th Edition

**Synonyms:**HsT1196;SCC;SCCA-1;SCCA-PD;SCCA1;T4-A

**Description:**SERPINB3, also known as SCCA-1, belongs to the serpin family. Serpins are a group of proteins with similar structures that were first identified as a set of proteins able to inhibit proteases. The acronym serpin was originally coined because many serpins inhibit chymotrypsin-like serine proteases. SERPINB3 is expressed in some hepatocellular carcinoma (at protein level). Its expression is closely related to cellular differentiation in both normal and malignant squamous cells. It seems to also be secreted in plasma by cancerous cells but at a low level. SERPINB3 significantly attenuates apoptosis by contrasting cytochrome c release from the mitochondria and by antichemotactic effect for NK cells. It may act as a protease inhibitor to modulate the host immune response against tumor cells and may be involved in the malignant behavior of squamous cell carcinoma cells.

**Form:**PBS

**Molecular Weight:**46.7 kDa

**Sequences:**Asn 2-Pro 390

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.