

Recombinant Mouse Cathepsin H Protein (His Tag)

Cat.NO.: TP07342

3th Edition

Synonyms:AL022844 Protein, Mouse

Description:Proinflammatory cytokines decreased the expression of CTSH in human islets and primary rat beta-cells, and overexpression of CTSH protected insulin-secreting cells against cytokine-induced apoptosis.

Form:PBS

Molecular Weight:36.6 kDa

Sequences:Ala21-Val333

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.