

Recombinant Rat UNC5B / UNC5H2 Protein (Fc tag)

Cat.NO.: TP07407

3th Edition

Synonyms:UNC5B

Description:The netrin-1 receptor, UNC-5 Homology B, or UNC5B plays vital roles in angiogenesis, inflammation, embryonic development and carcinogenesis. Overexpression of UNC5B human colon epithelial cells suppressed dextran sodium sulfate, or DSS-induced apoptosis and caspase-3 activity. In addition, is a potential anti-neoplastic target in bladder cancer progression and inflammatory arthritis.

Form:PBS

Molecular Weight:65.7 kDa

Sequences:Met1-Asp373

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.