

Recombinant Human IL17BR / IL17RB / IL-17 Receptor B Protein (His tag)

Cat.NO.: TP07442

3th Edition

Synonyms:CRL4;EVI27;IL17BR;IL17RH1

Description:MTSS1 (Metastasis suppressor 1), also known as Missing in metastasis (MIM), is a tissue-specific regulator of plasma membrane dynamics. MTSS1 is well described for its function as a metastasis suppressor gene and is expressed in a variety of tissues. MTSS1 might be involved in shaping neuronal membranes in vivo. MTSS1 deforms phosphoinositide-rich membranes through its I-BAR domain and interacts with actin monomers through its WH2 domain. MTSS1/MIM was first identified as a metastasis suppressor missing in metastatic bladder carcinoma cell lines. MTSS1 is a prognostic indicator of disease-free survival in breast cancer patients and demonstrates the ability to play a role in governing the metastatic nature of breast cancer cells. MTSS1 may serve as a useful biomarker for the prediction of outcome of gastric cancer. The down-regulation of MTSS1 that may be caused by DNA methylation was also observed in many other types of cancer. Recent work proposed that MIM also potentiates Sonic hedgehog (Shh)-induced gene expression. MTSS1 as a multiple functional molecular player and has an important role in development, carcinogenesis and metastasis.

Form:PBS

Molecular Weight:31.5 kDa

Sequences:Met 1-Gly289

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.