



Recombinant Rat Sclerostin / SOST Protein (His tag)

Cat.NO.: TP07495

3th Edition

Synonyms:SOST

Description:Sclerostin, the protein product of the SOST gene, is a potent inhibitor of bone formation. Sclerostin protein is widely expressed at low levels with highest levels in bone, cartilage, kidney, liver, bone marrow and primary osteoblasts differentiated for 21 days, and was originally identified as an important regulator of bone remodeling, homeostasis, and links bone resorption and bone apposition. Recent studies have revealed that Sclerostin protein inhibits the bone growth probably by binding to the extracellular domain of the Wnt coreceptors LRP5 and LRP6 and disrupting Wnt-induced Frizzled-LRP complex formation.

Form:PBS

Molecular Weight:23 kDa

Sequences:Met 1-Tyr 213

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.