

Recombinant Mouse KLRB1F / NKR-P1F Protein (Fc tag)

Cat.NO.: TP07520

3th Edition

Synonyms:A630024B12Rik;MGC130065;MGC130066;Nkr-p1f;Nkrp1-f;Nkrp1f

Description:

Form:PBS

Molecular Weight:46 kDa

Sequences:Gln 67-Val 217

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.