

Recombinant Mouse TCN2 / Transcobalamin-II Protein (His tag)

Cat.NO.: TP07599

3th Edition

Synonyms:AW208754;Tcn-2

Description:Transcobalamin II, also known as TCN2 and TC II, is a plasma protein that binds cobalamin (Cbl; vitamin B12) as it is absorbed in the terminal ileum and distributes to tissues. The circulating transcobalamin II-cobalamin complex binds to receptors on the plasma membrane of tissue cells and is then internalized by receptor-mediated endocytosis. Transcobalamin II is a non-glycosylated secretory protein of molecular mass 43 kDa. Its plasma membrane receptor (TC II-R) is a heavily glycosylated protein with a monomeric molecular mass of 62 kDa. Human TCN2 gene is composed of nine exons and eight introns spanning approximately 20 kb with multiple potential transcription start sites. A number of genetic abnormalities are characterized either by a failure to express TCN2 or by synthesis of an abnormal protein. The TCN2 deficiency results in cellular cobalamin deficiency, an early onset of megaloblastic anaemia, and neurological abnormalities.

Form:PBS

Molecular Weight:47 kDa

Sequences:Met 1-Trp 430

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.