

**Recombinant Mouse Notch-1 / NOTCH1 Protein (His tag)****Cat.NO.: TP07606**

3th Edition

**Synonyms:**9930111A19Rik;lin-12;Mis6;N1;Tan1

**Description:**NOTCH1 is one of the four mammalian Notch receptors, which is involved in the Notch signaling pathway. Specifically, NOTCH1 promotes the proliferation of myogenic precursor cells, and the NICD domain of NOTCH1 can impair regeneration of skeletal muscles. NOTCH1 is a prevalent signaling pathway in T cell acute lymphoblastic leukemia (T-ALL). The NOTCH signaling pathway is a conserved signaling cascade that regulates many aspects of development and homeostasis in multiple organ systems. The proto-oncogene NOTCH1 is frequently mutated in around 1% of patients with chronic lymphocytic leukemia (CLL). NOTCH1 mutations in oral squamous cell carcinoma (OSCC) frequently occur near the ligand-binding region. These mutations change the domain structure of this protein and affect the ligand binding activity. When NOTCH1 is activated by ligand binding, NOTCH1 intracellular domain (NICD) is cleaved from the cell membrane.

**Form:**PBS**Molecular Weight:**55 kDa**Sequences:**Met 1-Gln 526**Purity:**> 95% by HPLC**Concentration:****Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.