

---

**Recombinant Human CD16b / FCGR3B Protein (His tag)****Cat.NO.: TP07616**

3th Edition

**Synonyms:**CD16;CD16b;Fc gamma RIIIb;FCG3;FCGR3;FCR-10;FCRIII;FCRIIIb

**Description:**The asialoglycoprotein receptor (ASGPR), an endocytotic cell surface receptor expressed by hepatocytes, is triggered by triantennary binding to galactose residues of macromolecules such as asialoorosomuroid (ASOR). ASGPR belongs to the long-form subfamily of the C-type/Ca<sup>2+</sup> dependent lectin family. It is a complex of two noncovalently-linked and highly homologous subunits, a major 42 kDa glycoprotein ASGPR1(MHL-1) and a minor 51 kDa glycoprotein ASGR2 (MHL-2). ASGPR1 is synthesized as a type II transmembrane protein that contains a cytosolic N-terminal domain, a single transmembrane segment, and an extracellular domain which contains two important structural regions. The first is a stalk domain that contributes to noncovalent oligomerization, and the second is a Ca<sup>2+</sup>-dependent carbohydrate binding domain at the very C-terminus that is unusually stabilized by three ions. The research regarded that ASGPR1 could be targeted for anti-hepatitis B virus (HBV) drug development.

**Form:**PBS**Molecular Weight:**22.2 kDa**Sequences:**Met 1-Ser 200**Purity:**> 95% by HPLC**Concentration:****Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.