

---

**Recombinant Mouse S100A5 Protein (His tag)**

**Cat.NO.: TP07638**

3th Edition

**Synonyms:**S100D9

**Description:**S100 protein is a family of low molecular weight protein found in vertebrates characterized by two EF-hand calcium-binding motifs. There are at least 21 different S100 proteins, and the name is derived from the fact that the protein is 100% soluble in ammonium sulfate at neutral pH. Most S100 proteins are disulfide-linked homodimer, and is normally present in cells derived from the neural crest, chondrocytes, macrophages, dendritic cells, etc. S100 proteins have been implicated in a variety of intracellular and extracellular functions. They are involved in regulation of protein phosphorylation, transcription factors, the dynamics of cytoskeleton constituents, enzyme activities, cell growth and differentiation, and the inflammatory response. Protein S100-A5, also known as Protein S-100D, S100 calcium-binding protein A5, S100A5 and S100D, is a member of the S100 family which contains two EF-hand domains. S100A5 is also a novel member of the EF-hand superfamily of calcium-binding proteins that is poorly characterized at the protein level. It is expressed in very restricted regions of the adult brain. From birth onwards, S100A5 remained a neuronal-specific protein, only located in a subpopulation of neurons in the spiral ganglion.

**Form:**PBS

**Molecular Weight:**12.3 kDa

**Sequences:**Met 1-Lys 93

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.