

**Recombinant Human Transferrin Receptor / TFRC / CD71 Protein (His tag)**

**Cat.NO.: TP07639**

3th Edition

**Synonyms:**CD71;p90;T9;TFR;TFR1;TR;TRFR

**Description:**Transferrin receptor protein 1, also known as transferrin receptor, Trfr, p90, CD71 and TFRC, is a single-pass type II membrane protein which belongs to the peptidase M28 family and M28B subfamily. TFRC / CD71 is a membrane-bound protein expressed in larger amounts in proliferating. The specific expression of TFRC can represent a diagnostic tool or a therapeutic target in solid tumours expressing this antigen. Transferrin receptor is necessary for development of erythrocytes and the nervous system. TFRC / CD71 is regulated by cellular iron levels through binding of the iron regulatory proteins, IRP1 and IRP2, to iron-responsive elements in the 3'-UTR. Up-regulated upon mitogenic stimulation. TFRC / CD71 represents a marker of malignant transformation in the pancreas that could be applied as potential diagnostic and therapeutic target.

**Form:**PBS

**Molecular Weight:**77.4 kDa

**Sequences:**Cys 89-Phe 760

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.