

Mouse CHST5 / GlcNAc6ST-3 Protein (Fc Tag)

Cat.NO.: TP07662

3th Edition

Synonyms:AI173964;Gn6st-3;GST-4;I-GlcNAc-6-ST

Description:

Form:PBS

Molecular Weight:70.1 kDa

Sequences:Ser27-Ser395

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.