

**Mouse CHST5 / GlcNAc6ST-3 Protein (Fc Tag)**

**Cat.NO.: TP07662**

---

3th Edition

**Synonyms:** AI173964; Gn6st-3; GST-4; I-GlcNAc-6-ST

**Description:**

**Form:** PBS

**Molecular Weight:** 70.1 kDa

**Sequences:** Ser27-Ser395

**Purity:** > 95% by HPLC

**Concentration:**

**Endotoxin Level:** <1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:** Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.