

**Recombinant Mouse IL17BR / IL17RB / IL-17 Receptor B Protein (His tag)**

**Cat.NO.: TP07715**

3th Edition

**Synonyms:**Evi27;IL-17ER;IL-17Rh1;Il17br;IL17RH1

**Description:**

**Form:**PBS

**Molecular Weight:**31.3 kDa

**Sequences:**Met1-Gly286

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.