

Recombinant Mouse IL17BR / IL17RB / IL-17 Receptor B Protein (His tag)

Cat.NO.: TP07715

3th Edition

Synonyms:Evi27;IL-17ER;IL-17Rh1;Il17br;IL17RH1

Description:

Form:PBS

Molecular Weight:31.3 kDa

Sequences:Met1-Gly286

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.