

Recombinant Mouse EGFR / HER1 / ErbB1 Protein (Fc tag)**Cat.NO.: TP07758**

3th Edition

Synonyms:9030024J15Rik;A1552599;ErbB;Errb1;Errp;wa-2;wa2;Wa5

Description:As a member of the epidermal growth factor receptor (EGFR) family, EGFR protein is type I transmembrane glycoprotein that binds a subset of EGF family ligands including EGF, amphiregulin, TGF- β , betacellulin, etc. EGFR protein plays a crucial role in signaling pathway in the regulation of cell proliferation, survival and differentiation. Binding of a ligand induces EGFR protein homo- or heterodimerization, the subsequent tyrosine autophosphorylation and initiates various down stream pathways (MAPK, PI3K/PKB and STAT). In addition, EGFR signaling also has been shown to exert action on carcinogenesis and disease progression, and thus EGFR protein is proposed as a target for cancer therapy currently. Immune Checkpoint Immunotherapy Cancer Immunotherapy Targeted Therapy

Form:PBS**Molecular Weight:**96 kDa**Sequences:**Met 1-Ser 647**Purity:**> 95% by HPLC**Concentration:****Endotoxin Level:**<1.0 EU per 1 μ g of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.