

Recombinant Mouse SerpinA3c Protein (His tag)

Cat.NO.: TP07760

3th Edition

Synonyms:1A1;Kalbp;Klkbp;MGC144822;Serpina3c;spi2

Description:

Form:PBS

Molecular Weight:46 kDa

Sequences:Met 1-Ala 417

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.