

**Recombinant Human ITGA5 & ITGB1 Heterodimer Protein**

**Cat.NO.: TP07761**

3th Edition

**Synonyms:**CD49e;FNRA;VLA-5;VLA5A

**Description:**

**Form:**PBS

**Molecular Weight:**192 (109 + 83) kDa

**Sequences:**Met 1-Tyr 995

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.