



Recombinant Human ITGA5 & ITGB1 Heterodimer Protein

Cat.NO.: TP07761

3th Edition

Synonyms:CD49e;FNRA;VLA-5;VLA5A

Description:

Form:PBS

Molecular Weight:192 (109 + 83) kDa

Sequences:Met 1-Tyr 995

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.