

**Recombinant Mouse CES3 / Carboxylesterase-3 / CES1D Protein (His tag)**

**Cat.NO.: TP07765**

3th Edition

**Synonyms:**Ces1;Ces1d;CES3;TGH

**Description:**

**Form:**PBS

**Molecular Weight:**60.8 kDa

**Sequences:**Met1-Glu561

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.