

Recombinant Human ITGAX & ITGB2 Heterodimer Protein

Cat.NO.: TP07870

3th Edition

Synonyms:CD11C;SLEB6

Description:

Form:PBS

Molecular Weight:205 (125+80) kDa

Sequences:Met 1-Pro 1107

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.