

**Recombinant Mouse ALK-1 Protein (His &Fc Tag)**

**Cat.NO.: TP07877**

3th Edition

**Synonyms:**Acvrlk1 Protein, Mouse;Al115505 Protein, Mouse;Al427544 Protein, Mouse;Alk1 Protein, Mouse

**Description:**Activin A receptor, type II-like 1 (ACVRL1), also known as ALK-1 (activin receptor-like kinase 1), is an endothelial-specific type I receptor of the TGF-beta (transforming growth factor beta) receptor family of ligands. On ligand binding, a heteromeric receptor complex forms consisting of two type II and two type I transmembrane serine/threonine kinases. ACVRL1 protein is expressed in certain blood vessels of kidney, spleen, heart and intestine, serving as an important role during vascular development. Mutations in ACVRL1 gene are associated with hemorrhagic telangiectasia type 2, also known as Rendu-Osler-Weber syndrome 2 and vascular disease.

**Form:**PBS

**Molecular Weight:**39 kDa

**Sequences:**Met 1-Pro 119

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.