

Recombinant Mouse PILR-alpha / PILRA Protein (His tag)

Cat.NO.: TP07918

3th Edition

Synonyms:AV021745;FDF03

Description:

Form:PBS

Molecular Weight:20.2 kDa

Sequences:Met 1-Val 197

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.