

Recombinant Human Purine nucleoside phosphorylase / PNP Protein (His tag)

Cat.NO.: TP07929

3th Edition

Synonyms:NP;PRO1837;PUNP

Description:RGMa, also known as RGM domain family, member A, belongs to the RGM (repulsive guidance molecule) family whose members are membrane-associated glycoprotein. RGMa is a glycosylphosphatidylinositol-anchored glycoprotein that functions as an axon guidance protein in the developing and adult central nervous system. It helps guide Retinal Ganglion Cell (RGC) axons to the tectum in the midbrain. RGMa has been implicated to play an important role in the developing brain and in the scar tissue that forms after a brain injury. This protein may also function as a tumor suppressor in some cancers.

Form:PBS

Molecular Weight:33.5 kDa

Sequences:Met 1-Ser 289

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.