

**Recombinant Human Isocitrate dehydrogenase Protein (His Tag)**

**Cat.NO.: TP08018**

3th Edition

**Synonyms:**HEL-216;HEL-S-26;IDCD;IDH;IDP;IDPC;PICD

**Description:**

**Form:**PBS

**Molecular Weight:**48 kDa

**Sequences:**Met 1-Leu 414

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.