

Recombinant Mouse CD16-2 / FCGR4 Protein (His tag)

Cat.NO.: TP08025

3th Edition

Synonyms:4833442P21Rik;CD16-2;FcgammaRIV;Fcgr3a;FcgRIV;Fcrl3

Description:

Form:PBS

Molecular Weight:22.4 kDa

Sequences:Met 1-Gln 203

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.