

**Recombinant Mouse Neurexophilin-1 / NXPH1 Protein (His tag)**

**Cat.NO.: TP08043**

3th Edition

**Synonyms:**C130005L03Rik

**Description:**Neurexophilin-1, or NXPH1 is a secreted glycoprotein, which belongs to the Neurexophilin family. The Neurexophilin family contain at least four genes and resembles a neuropeptide, suggesting a function as an endogenous ligand for alpha-neurexins. The mammalian brains contain four genes for neurexophilins the products of which share a common structure composed of five domains: an N-terminal signal peptide, a variable N-terminal domain, a highly conserved central domain that is N-glycosylated, a short linker region, and a conserved C-terminal domain that is cysteine-rich. Neurexophilin-1 constitutes a secreted cysteine-rich glycoprotein, forms a very tight complex with alpha neurexins, a group of proteins that promote adhesion between dendrites and axons. Neurexophilins 1 and 3 but not 4 (neurexophilin 2 is not expressed in rodents) bind to a single individual LNS domain, the second overall LNS domain in all three alpha-neurexins.

**Form:**PBS

**Molecular Weight:**30.1 kDa

**Sequences:**Ala22-Gly271

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.