

---

**Recombinant Mouse OX40 Protein (Fc Tag)****Cat.NO.: TP08102**

3th Edition

**Synonyms:**ACT35 Protein, Mouse;CD134 Protein, Mouse;Ly-70 Protein, Mouse;Ox40 Protein, Mouse;Txgp1 Protein, Mouse;TXGP1L Protein, Mouse

**Description:**OX40 (CD134) and its binding partner, OX40L (CD252), are members of the tumor necrosis factor receptor/tumor necrosis factor superfamily, is known to break an existing state of tolerance in malignancies, leading to a reactivation of antitumor immunity. The interaction between OX40 and OX40L plays an important role in antigen-specific T-cell expansion and survival. OX40 and OX40L also regulate cytokine production from T cells, antigen-presenting cells, natural killer cells, and natural killer T cells, and modulate cytokine receptor signaling. In line with these important modulatory functions, OX40-OX40L interactions have been found to play a central role in the development of multiple inflammatory and autoimmune diseases, making them attractive candidates for intervention in the clinic. Conversely, stimulating OX40 has shown it to be a candidate for therapeutic immunization strategies for cancer and infectious disease.

**Form:**PBS

**Molecular Weight:**48.3 kDa

**Sequences:**Met 1-Pro 211

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.