
Recombinant Human Serglycin / SRGN Protein (His & Myc tag)**Cat.NO.: TP08110**

3th Edition

Synonyms:PPG;PRG;PRG1

Description:SRGN is known as a hematopoietic cell granule proteoglycan. Proteoglycans stored in the secretory granules of various hematopoietic cells has a protease-resistant peptide core, and is vital for neutralizing hydrolytic enzymes. SRGN is associated with the macromolecular complex of granzymes and perforin, which may serve as a mediator of granule-mediated apoptosis. It is required for storage of some proteases in both connective tissue and mucosal mast cells and for storage of granzyme B in T-lymphocytes. SRGN also plays a role in localizing neutrophil elastase in azurophil granules of neutrophils.

Form:PBS**Molecular Weight:**17.3 kDa**Sequences:**Met 1-Leu158**Purity:**> 95% by HPLC**Concentration:****Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.