

Recombinant Mouse LAMP2 / CD107b Protein (His tag)

Cat.NO.: TP08130

3th Edition

Synonyms:CD107b;Lamp-2;Lamp-2a;Lamp-2b;Lamp-2c;LampII;LGP-B;Mac3

Description:LAMP2 (Lysosomal-associated membrane protein 2), also known as CD107b (Cluster of Differentiation 107b), is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. In human, LAMP2, the causative gene of Danon disease, located on chromosome Xq24, encodes the lysosome-associated membrane protein-2 (LAMP-2). LAMP-2 deficiency, or Danon disease, is a rare X-linked lysosomal disease characterized by cardiomyopathy, vacuolar myopathy, and mental retardation. LAMP2 cardiomyopathy is an X-linked and highly progressive myocardial storage disorder associated with diminished survival, which clinically resembles sarcomeric hypertrophic cardiomyopathy.

Form:PBS

Molecular Weight:40.6 kDa

Sequences:Leu 26-Asn 379

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.