
Recombinant Rat c-MPL / CD110 / TPOR Protein (His tag)**Cat.NO.: TP08132**

3th Edition

Synonyms:MPL

Description:CD110, also known as c-MPL, is a 635 amino acid transmembrane domain, with two extracellular cytokine receptor domains and two intracellular cytokine receptor box motifs. It is expressed at a low level in a large number of cells of hematopoietic origin. C-MPL is homologous with members of the hematopoietic receptor superfamily. Presence of anti-sense oligodeoxynucleotides of c-mpl inhibited megakaryocyte colony formation. Thrombopoietin is the ligand for c-mpl. It was shown to be the major regulator of megakaryocytopoiesis and platelet formation. Defects in c-MPL are a cause of congenital amegakaryocytic thrombocytopeniawhich is a disease characterized by isolated thrombocytopenia and megakaryocytopenia with no physical anomalies. Defects in c-MPL also cause thrombocythemia type 2 and myelofibrosis with myeloid metaplasia.

Form:PBS**Molecular Weight:**55.1 kDa**Sequences:**Met1-Ala500**Purity:**> 95% by HPLC**Concentration:****Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.