
Recombinant Human SerpinA3 / AACT Protein (His tag)**Cat.NO.: TP08143**

3th Edition

Synonyms:AACT;ACT;GIG24;GIG25;MGC88254;SERPINA3

Description:SerpinA3, also known as Alpha 1-antichymotrypsin (AACT), is a plasma alpha globulin glycoprotein, and is a member of serpin superfamily of the serine protease inhibitors consisting of at least 35 members. SerpinA3 has been demonstrated to inhibit the activity of certain serine proteases, such as cathepsin G found in neutrophils, and chymases present in mast cells, by inducing a major conformational rearrangement, and thus protects some tissues from damage caused by proteolytic enzymes. This enzyme is produced primarily in the liver, and is identified as an acute-phase inflammatory protein. SerpinA3 deficiency has been associated with liver disease, and mutations of this gene have been observed in patients with Parkinson disease and chronic obstructive pulmonary disease. In addition, ACT gene polymorphism has been implicated with Alzheimer's disease (AD), cerebral amyloid angiopathy (CAA), as well as stroke, since SerpinA3 is a major constituent of the plaques in AD and an inhibitor of amyloid beta peptide degradation.

Form:PBS**Molecular Weight:**46.5 kDa**Sequences:**Met 1-Ala 423**Purity:**> 95% by HPLC**Concentration:****Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.