

Recombinant Human IL17BR / IL17RB / IL-17 Receptor B Protein (Fc tag)

Cat.NO.: TP08160

3th Edition

Synonyms:CRL4;EVI27;IL17BR;IL17RH1

Description:R-spondin-2, also known as RSPO2, synergizes with Wnt to activate beta-catenin. RSPO2 is secreted proteins that regulate beta-catenin signaling. Activator of the beta-catenin signaling cascade leads to TCF-dependent gene activation. Action both in the canonical Wnt / beta- catenin-dependent pathway, possibly via a direct interaction with Wnt proteins, and in a Wnt-independent beta catenin pathway through a receptor signaling pathway that may not use frizzled / LRP receptors. Probably also acts as a ligand for frizzled and LRP receptors. The encoding gene Rspo2 was identified as a novel common integration site for the mouse mammary tumor virus in viral induced mouse mammary tumors. Rspo2 and Rspo2 / Wnt1 tumors contained many spindle cells, consistent with an epithelial-mesenchymal transformation phenotype. When Rspo2 and Rspo2 / Wnt1 tumor cells were transferred into naive mice, they exhibited greater metastatic activity than cells derived from Wnt1 tumors.

Form:PBS

Molecular Weight:56.9 kDa

Sequences:Met 1-Gly289

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.