

**Recombinant Mouse CD16-2 / FCGR4 Protein (His & AVI tag)**

**Cat.NO.: TP08278**

3th Edition

**Synonyms:**4833442P21Rik;CD16-2;FcgammaRIV;Fcgr3a;FcgRIV;Fcrl3

**Description:**

**Form:**PBS

**Molecular Weight:**24.2 kDa

**Sequences:**Met 1-Gln 203

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.