

**Recombinant Human CRELD1 Protein (His tag)**

**Cat.NO.: TP08281**

3th Edition

**Synonyms:**AVSD2;CIRPIN

**Description:**CRELD1 is a transmembrane glycoprotein. Epidermal growth factor(EGF)like domain exists in CRELD1. EGF-like repeats are a class of cysteine-rich domains that mediate interactions between proteins of diverse function. EGF domains are found in proteins that are either completely secreted or have transmembrane regions that tether the protein to the cell surface. CRELD1 contains a 333 amino acid (aa) extracellular domain (ECD), two tandem transmembrane segments, and a second ECD of 15 aa. Defects in CRELD1 may cause susceptibility to atrioventricular septal defect type 2 which results in a persistent common atrioventricular canal.

**Form:**PBS

**Molecular Weight:**37.4 kDa

**Sequences:**Met 1-Glu362

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.