

Recombinant Rat FCGRT & B2M Heterodimer Protein

Cat.NO.: TP08516

3th Edition

Synonyms:FCGRT FCRN & B2M

Description:

Form:PBS

Molecular Weight:44 kDa

Sequences:Met 1-Ser 298

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.