



Recombinant Human AMPK (G1/B1/A2) Heterotrimer Protein

Cat.NO.: TP08628

3th Edition

Synonyms:AMPK;AMPKa1

Description:

Form:PBS

Molecular Weight:158 (38+30+90) kDa

Sequences:Met 1-Pro 331

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.