

**Recombinant Human AMPK (G1/B2/A2) Heterotrimer Protein**

**Cat.NO.: TP08630**

3th Edition

**Synonyms:**AMPK;AMPKa1

**Description:**

**Form:**PBS

**Molecular Weight:**158 (38+30+90) kDa

**Sequences:**Met 1-Pro 331

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.