

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ACTB Anti-Human Antibody (Clone AT4G6)

单克隆抗体

产品货号: MA01009

第三版

免疫原: Recombinant human ACTB (1-375aa) purified from E. coli.

亚型: Mouse IgG2a heavy chain and kappa light chain

克隆: Anti-human ACTB mAb, clone AT4G6, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ACTB protein.

描述: Beta-actin (human gene and protein symbol ACTB/ACTB) is one of six different actin isoforms which have been identified in humans. This is one of the two nonmuscle cytoskeletal actins. Actins are highly conserved proteins that are involved in cell motility, structure and integrity. Alpha actins are a major constituent of the contractile apparatus.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF